

## HIGH PRECISION TABLE TOP SCALE



### SALIENT FEATURES

- Fast & Stable Display
- Single Touch Key Calibration with external weights
- Easy to Read LED Display
- Multifunction Weighing units
- Auto Zero Tracking
- Overload Protection Design
- Stainless Steel Weighing Pan
- Adjustable Feet
- High Resolution
- Quick & Accurate Result

## SPECIFICATIONS

MODEL	SPT 3C	SPT 6C	SPT 15C	SPT 30C
Capacity	3 kg	6 kg	15 kg	30 kg
Readability	0.1 g	0.2 g	0.5 g	1 g
Repeatability	0.1 g	0.2 g	0.5 g	1 g
Pan Size	255 x 190 mm			
Response Time	3 – 4 seconds			
Calibration	Automatic Calibration with External Weight			
Tare Range	100% of F.S.			
Operating Temperature	10 c to 40 c			
Dimensions				
Power Supply	220V / 50-60 Hz			
Net / Gross Weight	Approx 6 kg			