



LAB BALANCE



SALIENT FEATURES

- 2,00,000 Internal Resolutions
- Auto Calibration with external weights
- Easy to Read Big LCD Display with blue backlight
- Multifunction Weighing units
- Auto Zero Tracking
- Aluminium casting body
- Overload Protection Design
- Stainless Steel Weighing Pan
- Adjustable Feet
- High Resolution
- Quick & Accurate Result

SPECIFICATIONS

MODEL	SPT 200 CK	SPT 300 CK
Capacity	200 g	300 g
Readability	0.001 g	0.001 g
Repeatability	0.001 g	0.001 g
Linearity ±	0.002 g	
Pan Size	85 mm dia	
Response Time	3 – 4 seconds	
Calibration	Automatic Calibration with External Weight	
Tare Range	100% of F.S.	
Operating Temperature	10 c to 40 c	
Dimensions		
Power Supply	AC adapter 220V / 50-60 Hz	
Net / Gross Weight	Approx 6 kg	

Prime Technologies

935, Sector 15A, Faridabad – 121 007, Haryana, INDIA.

Telephone : +91 129 2222698, 4015935

Mobile : +91 93501 85341 (WhatsApp)

E-mail : primetechnoindia@gmail.com | mahavirweightronix@gmail.com

Web : www.primetechnoindia.com