




## LAB BALANCE (JC SERIES)

### SALIENT FEATURES



- 2,00,000 Internal Resolutions
- Auto Calibration with external weights
- Easy to Read Big LCD Display with white backlight
- Multifunction Weighing units
- Auto Zero Tracking
- Low Battery Indication
- Overload Protection Design
- Stainless Steel Weighing Pan
- Adjustable Feet
- High Resolution
- Quick & Accurate Result
- In built Battery backup

### SPECIFICATIONS

MODEL	SPT 120 JC	SPT 220 JC	SPT 320 JC
Capacity	120 g	220 g	320 g
Readability	0.001 g	0.001 g	0.001 g
Repeatability	0.001 g	0.001 g	0.001 g
Linearity ±	0.002 g		
Pan Size	85 mm dia		
Response Time	3 – 4 seconds		
Calibration	Automatic Calibration with External Weight		
Tare Range	100% of F.S.		
Operating Temperature	10 c to 40 c		
Dimensions	280 x 188 x 290 mm		
Power Supply	AC adapter 220V / 50-60 Hz		
Net / Gross Weight	Approx 4.1 kg		

## Prime Technologies

935, Sector 15A, Faridabad – 121 007, Haryana, INDIA.

Telephone : +91 129 2222698, 4015935

Mobile : +91 93501 85341 (WhatsApp)

E-mail : [primetechnoindia@gmail.com](mailto:primetechnoindia@gmail.com) | [mahavirweightronix@gmail.com](mailto:mahavirweightronix@gmail.com)

Web : [www.primetechnoindia.com](http://www.primetechnoindia.com)