




MONOLITHIC MAGNETIC LAB BALANCE



SALIENT FEATURES

- Based on Monolithic Magnetic Force Technology for high Weighing accuracy & reliability for daily lab routines
- Single Touch Key Calibration with external weights
- Easy to Read high contrast & Big LCD Display with backlight
- Multifunction Weighing units
- Auto Zero Tracking
- Overload Protection Design
- Stainless Steel Weighing Pan
- Adjustable Feet
- High Resolution
- Glass Draft Shield
- Parts Counting Mode

SPECIFICATIONS

MODEL	SPT 120 ME	SPT 220 ME
Capacity	120 g	220 g
Readability	0.1 mg	0.1 mg
Repeatability	0.2 mg	0.2 mg
Linearity ±	0.3 mg	
Pan Size	SS weighing pan with removable dust & spill rings ~ 80 mm dia	
Response Time	< 5 seconds	
Calibration	Automatic Single Touch Key Calibration with External Weight	
Tare Range	100% of F.S.	
Operating Temperature	10 c to 40 c	
Dimensions		
Power Supply	AC adapter 220V / 50-60 Hz	
Net / Gross Weight	Approx 6 kg	

Prime Technologies

935, Sector 15A, Faridabad – 121 007, Haryana, INDIA.

Telephone : +91 129 2222698, 4015935

Mobile : +91 93501 85341 (WhatsApp)

E-mail : primetechnoindia@gmail.com | mahavirweightronix@gmail.com

Web : www.primetechnoindia.com